




Culprit Coronary Artery Algorithm (CCAA)

Culprit Artery Algorithm: Because time matters



SCHILLER

The Art of Diagnostics

SCHILLER CCAA

Software that helps correctly identify the site of obstruction in the coronary artery

SCHILLER CCAA points the way to the optimal treatment of acute chest pain, benefiting the patient and ensuring best use of hospital facilities.

Our objective

- to use the ECG to determine the size of the cardiac area at risk by localizing the occlusion site in the coronary artery
- provide clinical data to shorten the time interval between onset of chest pain and restoration of myocardial blood flow

The need

- shorten the time interval between onset of pain and medical attention
- accurate and quick decision-making using the standard 12-lead ECG

Our new approach

- Uses the ST-segment deviation direction of 12 ECG leads to indicate the site of occlusion in the culprit artery
- That information is combined with other ECG findings and a limited number of questions
- A software algorithm including a decision tree provides suggestions about the preferred therapy


In the US approximately 40 % of all deaths are caused by cardiovascular disease.

A decision about the correct treatment of the patient with acute chest pain can be made, while outside hospital by the ambulance personnel, in the hospital emergency room, or at home by the general practitioner.

Crucial time is lost in risk stratifying patients with acute chest pain

In many countries the number of hospitals where a percutaneous coronary intervention (PCI) can be performed is limited. This stresses the necessity of rapidly determining which patient should be brought to such a hospital.

The Challenge



	I	II	III	AVR	AVL	AVF
ST Values(uV) :	-28	-100	-73	65	22	-8
ST Values(mm) :	0.3	-1.0	0.7	0.7	0.2	0.1
ST Score (mm) :	6.8					
Score table :	I	II	III	AVR	AVL	AVF
Left main :	1	1	1	1	0	0
LAD Prox :	1	1	1	1	1	1
LAD Dist :	1	0	1	0	1	1
RCA Prox :	1	0	0	0	1	1
RCA Dist :	1	0	0	0	1	1
Cx :	1	0	1	0	1	1
3V/LM Nar. :	1	1	1	1	0	0
Previous MI, CA Bypass or Stenting :	YES					
Time interval onset chest pain and ECG :	1.5 h					
QRS width :	82 ms					
ST Score :	6.8 mm					
Occlusion Site:						
- LAD Prox						
Advice:						
- PCI Center						
- Consider thrombolytic therapy if PCI centre is further						

The optional SCHILLER CCAA software, developed by Professor Hein Wellens and implemented in our ECG machines, allows the identification of the site of obstruction in the coronary artery.

This revolutionary algorithm points the way to the optimal treatment of acute chest pain, benefiting the patient, staff and hospital.

Effectively use the 12-lead ECG to optimize diagnosis and management of the patients with acute chest pain

12-lead ECG recording with lead V4R instead of lead V4. Lead V4R is a very helpful lead in infero-posterior myocardial infarction to distinguish between a coronary occlusion site proximal or distal in the right coronary artery, or in the circumflex coronary artery.

It has been shown that the ST-segment

SCHILLER CCAA

Helps ensure a correct treatment decision


Essential for correct 12-lead ECG interpretation and conformation of STEMI events

- Correct ECG lead positioning
- No artifacts
- Localization of the ischemic area with accurate determination of site and size of the coronary artery occlusion
- Identification of other ECG and non-ECG risk factors

Benefits of the SCHILLER CCAA algorithm

- Saving lives
- Shortens hospital stay and transport costs
- Cost-effective, quick pre-hospital decision-making regarding next treatment steps
- Appropriate selection and triage of patients requiring intra-coronary intervention / thrombolysis
- Reducing the duration and size of the myocardial infarction and damage to the heart muscle

Literature: Meissner A, Trappe H-J, de Boer M-J, Gorgels AP, Wellens HJ. «The value of the ECG for decision-making at first medical contact in the patient with acute chest pain». Neth Heart J. 2010;18:301-306



The closer the occlusion site to the origin of the coronary artery, the larger the size of the area at risk

	V1	V2	V3	V4	V5	V6	V4R
F	62	47	23	-44	-65	-63	-18
S	0.6	0.5	0.2	0.4	0.7	0.6	


	V1	V2	V3	V4	V5	V6	V4R	->	Sum
1	1	0	1	1	1	1	0	->	10
1	1	0	0	1	1	1	1	->	11
1	0	0	0	1	1	1	1	->	8
1	1	1	1	1	1	1	0	->	9
0	1	0	1	1	1	1	1	->	8
1	0	0	0	1	1	1	0	->	7
1	1	0	1	1	1	1	0	->	10

other away than 1.5 hours.


behavior of lead V4 can be derived by averaging the ST-segment deviation of leads V3 and V5. Correct placement of lead V4R as

indicated in the figure below. It is used in the Wellens algorithm for added efficacy. The algorithm indicates the location

of the different occlusion sites in the coronary artery as shown in the diagram below.



RCA: Right Coronary artery
RV branch: Right Ventricular
LM: Left Main
RCX: Ramus Circumflexus
LAD: Left Anterior Descending





SCHILLER, founded in 1974 by Alfred E. Schiller is a successful group with 31 subsidiaries and a global sales network. Today, SCHILLER is a world-leading manufacturer and supplier of devices for cardiopulmonary diagnostics, defibrillation and patient monitoring as well as software solutions for the medical industry.

Over the past 25+ years, Schiller India has established itself as a leader in the medical technology sector. With 450+ employees, a state-of-the-art ISO 13485-certified production centre in Puducherry, and multiple R&D centres across the globe, Schiller India makes advanced healthcare equipment accessible through a network of 100 sales and service dealers across more than 45 locations. Our product range includes Critical Care, Anaesthesia, Emergency Care, Cardiology, Respiratory Diagnostics, Radiology, and Robotics.

All registered trademarks acknowledged.

SCHILLER
The Art of Diagnostics

Swiss H.Q.: SCHILLER AG, Altgasse 68, P. O. Box 1052, CH -6341 Baar, Switzerland

Indian H.O.: Schiller Healthcare India Pvt. Ltd., Advance House, Makwana Road, Andheri (East), Mumbai-59. Tel.: + 91 - 22 - 61523333

Factory: No. 15/5 & 15/6, Vazhuthavur Road, Kurumbapet, Puducherry - 605 009, India

Regional Offices: Mumbai: 022-61523333 Email: sales.west@schillerindia.com Delhi: 011-41062067 / 09312432205 Email: sales.north@schillerindia.com Kolkata: 033-23593102 / 033-23593103 Email: sales.east@schillerindia.com Chennai: 044-28232648 / 044-28311021 Email: sales.south@schillerindia.com Ahmedabad: 079- 35337949 / 079- 35337950 Email: sales.west@schillerindia.com Puducherry: 0413- 220 2680 / 0413- 229 2940 Email: sales.south@schillerindia.com

CIN : U33110MH1997PTC111307

sales@schillerindia.com, help@schillerindia.com www.schillerindia.com 1800 209 8998

Follow us on : [f/SchillerHealthcareIndia](#) [i/schiller_healthcare_india](#) [in/schillerhealthcareindia](#) [s/schiller_healthcare_india](#) [@schillerindia1074](#)