

SpiroScout SP plus

PC Spirometer

Ultrasound Spirometry measurements with CS-104 software



SPIROSCOUT SP PLUS

SpiroScout SP plus:

- ❖ Accurate GANSHORN ultrasound technology
- ❖ Calibration-free
- ❖ No warm-up time and no moving parts

The SpiroScout SP plus incorporates an ultrasound flow sensor to measure the flow of air into and out of the patients lungs. The ultrasound technology eliminates problems associated with traditional methods of flow measurement: there are no moving parts, no screens to catch sputum and no disposables to calibrate.

Ultrasound flow measurement is independent of gas composition, pressure and humidity and it eliminates errors connected with these variables. The disposable ScoutTube lets the ultrasonic pulses travel between the transducers.

The disposable ScoutTube does not perform any measurement functions and therefore does not require calibration.

DISPOSABLE AND SUSTAINABLE

To ensure a sustainable product, SCHILLER compensates the CO₂ emitted during production, packaging, shipment, and disposal of the Scout-Tubes and the PFT filters.



01-20-162037
myclimate.org

UNIQUE ADVANTAGES OF GANSHORN ULTRASOUND TECHNOLOGY:


Advantages of Ultrasound Technology

- ❖ Sensor not in contact with sample
- ❖ Not influenced by humidity, barometric pressure, contamination
- ❖ Calibration-free
- ❖ No cleaning of sensors / filters, maintenance-free
- ❖ Simple and economic consumables
- ❖ Extremely high accuracy for low flows
- ❖ No downtime

Disadvantages of Traditional Technology*

- ❖ Sensor always in contact with sample
- ❖ Extremely sensitive to contamination (measuring errors)
- ❖ Regular calibration needed
- ❖ Need for cleaning
- ❖ Limited low-flow detection
- ❖ Downtime for maintenance

*turbine, pneumotach, heated wire



SCOUTTUBE

- ❖ Hygienic single-patient consumable for full contamination control
- ❖ Oval in shape for optimal lip seal and teeth positioning
- ❖ Conical to suit different mouth sizes

PFT FILTER

- ❖ Provides accurate, leakage-free and reproducible results
- ❖ Protects against cross-contamination and offers additional protection for operators and technicians
- ❖ Disinfectable adapter



About SCHILLER

SCHILLER A. G., a company based in Switzerland, was founded by Alfred Schiller in 1974. Over the last 50 years, it has evolved into a successful Global Group with 31 subsidiaries and dealers across 30 countries. Today, SCHILLER is a world-leading manufacturer and supplier of devices for cardiopulmonary diagnostics, defibrillation and patient monitoring as well as software solutions for the medical industry.

Schiller in India was founded as a joint venture in 1997, and since then has established itself as a leader in the medical technology sector. With 450+ employees, a state-of-the-art ISO 13485-certified production centre in Puducherry, multiple R&D centres, and 17 offices across India; Schiller makes advanced healthcare equipment accessible through a network of 100 sales and service dealers across more than 45 locations. Our product range includes Critical Care, Anaesthesia, Emergency Care, Cardiology, Respiratory Diagnostics, Radiology, and Robotics.

Scan here
to know more:



All registered trademarks acknowledged.

SCHILLER
The Art of Diagnostics

Swiss H.Q.: SCHILLER AG, Altgasse 68, P. O. Box 1052, CH -6341 Baar, Switzerland

Indian H.O.: Schiller Healthcare India Pvt. Ltd., Advance House, Makwana Road, Andheri (East), Mumbai-59. Tel : + 91 - 22 - 61523333

: No. 15/5 & 15/6, Vazhuthavur Road, Kurumbapet, Puducherry - 605 009, India

Regional Offices: Mumbai: 022-61523333 Email: sales.west@schillerindia.com Delhi: 011-41062067 / 09312432205 Email: sales.north@schillerindia.com Kolkata: 033-23593102 / 033-23593103 Email: sales.east@schillerindia.com Chennai: 044-28232648 / 044-28311021 Email: sales.south@schillerindia.com Ahmedabad: 079-35337949 / 079-35337950 Email: sales.west2@schillerindia.com Puducherry: 0413-2202680 / 0413-2292940 Email: sales.south@schillerindia.com

CIN : U33110MH1997PIC111307

sales@schillerindia.com, help@schillerindia.com www.schillerindia.com 1800 209 8998

Follow us on : /SchillerHealthcareIndia /schiller_healthcare_india /schillerhealthcareindia /schiller_healthcare_india @schillerindia1074