

TRUSCOPE[®] III

12.1" screen multi-utility monitor



SCHILLER
The Art of Diagnostics

TRUSCOPE[®] III

TRUSCOPE III monitor is a high-performance patient monitoring system designed to deliver accurate and real-time vital sign measurements, ensuring reliable care in critical environments.



Advanced Patient Monitoring Features

- ✦ TRUSCOPE III features a 12.1" high-resolution color TFT LCD display, ensuring clear and detailed visualizations.
- ✦ With 168 hours of tabular and graphic trend data storage, healthcare professionals can easily track patient progress over time.
- ✦ The standard configuration includes ECG, Respiration (RESP), Non-invasive Blood Pressure (NIBP), SpO₂ and Temperature monitoring.
- ✦ It is equipped with a pacemaker detection facility and offers two USB ports for enhanced connectivity and data transfer.



User Friendly Design

- ✦ The TRUSCOPE III monitor features an intuitive interface with customizable smart keys, simplifying complex tasks and streamlining workflows.
- ✦ This design boosts efficiency, minimizes errors, and helps healthcare professionals concentrate on patient care during high-pressure situations.

Comprehensive Network Monitoring

- ✦ TRUSCOPE III ensures seamless connectivity, offering both wired and Wi-Fi options to maintain a continuous and reliable link to the central nursing station.
- ✦ Enables healthcare professionals to monitor vital sign parameters in real-time, ensuring prompt and accurate patient care.
- ✦ The system's full compatibility with HIS / HL-7 standards, enhancing integration and ensuring efficient communication across all healthcare platforms.

Effortless Data Management and Transfer


- ❖ Our monitoring system is designed with user convenience in mind, offering straightforward data management solutions.
- ❖ The USB export feature simplifies the transfer of trend data for accessible critical patient information
- ❖ The plug-and-play interface supports a variety of peripheral devices, such as USB sticks and mice, making data retrieval and device operation both intuitive and efficient.
- ❖ The system supports data export directly to a PC, which can be viewed and analyzed using our dedicated PC Viewer Software, providing healthcare professionals with the tools they need to make informed decisions.




Customizable Screen Layouts for Optimal Viewing

- ❖ Understanding the diverse needs of healthcare environments, our system comes equipped with eight different screen layouts.
- ❖ Customizable views are designed to accommodate various angles and distances, ensuring that critical information is always visible and accessible, no matter where the healthcare provider is positioned.
- ❖ Enhances workflow efficiency and allows for more flexible monitoring solutions tailored to the specific needs of each user.


Standard



7- Lead



12- Lead



Big Numerics





About SCHILLER

SCHILLER A. G., a company based in Switzerland, was founded by Alfred Schiller in 1974. Over the last 50 years, it has evolved into a successful Global Group with 31 subsidiaries and dealers across 30 countries. Today, SCHILLER is a world-leading manufacturer and supplier of devices for cardiopulmonary diagnostics, defibrillation and patient monitoring as well as software solutions for the medical industry.

Schiller in India was founded as a joint venture in 1997, and since then has established itself as a leader in the medical technology sector. With 450+ employees, a state-of-the-art ISO 13485-certified production centre in Puducherry, multiple R&D centres, and 17 offices across India; Schiller makes advanced healthcare equipment accessible through a network of 100 sales and service dealers across more than 45 locations. Our product range includes Critical Care, Anaesthesia, Emergency Care, Cardiology, Respiratory Diagnostics, Radiology, and Robotics.



Scan here
to know more:

Part no.: B2.100005 09/2025 Ver. no.: 1.2

All registered trademarks acknowledged.

SCHILLER
The Art of Diagnostics

Swiss H.Q.: SCHILLER AG, Altgasse 68, P. O. Box 1052, CH -6341 Baar, Switzerland

Indian H.O.: Schiller Healthcare India Pvt. Ltd., Advance House, Makwana Road, Andheri (East), Mumbai-59. Tel : + 91 - 22 - 61523333

Factory: No. 15/5 & 15/6, Vazhuthavur Road, Kurumbapet, Puducherry - 605 009, India

Regional Offices: Mumbai: 022-61523333 Email: sales.west@schillerindia.com Delhi: 011-41062067 / 09312432205 Email: sales.north@schillerindia.com Kolkata: 033-23593102 / 033-23593103 Email: sales.east@schillerindia.com Chennai: 044-28232648 / 044-28311021 Email: sales.south@schillerindia.com Ahmedabad: 079- 35337949 / 079- 35337950 Email: sales.west@schillerindia.com Puducherry: 0413- 220 2680 / 0413- 229 2940 Email: sales.south@schillerindia.com

CIN : U33110MH1997PTC111307

sales@schillerindia.com, help@schillerindia.com www.schillerindia.com 1800 209 8998

Follow us on : /SchillerHealthcareIndia /schiller_healthcare_india /schillerhealthcareindia /schiller_healthcare_india @schillerindia1074