

MUSTANG

The MUSTANG is a highly durable and silent treadmill, with a spacious walking deck. Its smooth speed and inclination control makes it ideal for cardiopulmonary testing and cardiac rehabilitation.



ENHANCED USABILITY AND PERFORMANCE



High-Performance Motor: Powered by a 2.0 HP heavy-duty AC motor with a speed range of 0-22 Km/h and an incline range of 0-25% grade, enabling rigorous stress testing protocols.



Spacious Running Surface: A generous 1500 mm x 500 mm walking / running area ensures comfortable movement for users weighing up to 250 Kg.



Quiet and Smooth Performance: Equipped with advanced motor technology, it ensures low noise and a smooth, joint-friendly workout experience at any speed.



Safe & Compliant Design: Equipped with an emergency stop switch, Type B shock protection, and compliance with safety standards for electrical leakage and ingress protection.



Stable & Low-maintenance: It has a strong frame, is easy to maintain, and provides a smooth and comfortable run because of its smart design and powerful motor.



SCHILLER
The Art of Diagnostics



SCHILLER

The Art of Diagnostics

About SCHILLER

SCHILLER A. G., a company based in Switzerland, was founded by Alfred Schiller in 1974. Over the last 50 years, it has evolved into a successful Global Group with 31 subsidiaries and dealers across 30 countries. Today, SCHILLER is a world-leading manufacturer and supplier of devices for cardiopulmonary diagnostics, defibrillation and patient monitoring as well as software solutions for the medical industry.

Schiller in India was founded as a joint venture in 1997, and since then has established itself as a leader in the medical technology sector. With 450+ employees, a state-of-the-art ISO 13485-certified production centre in Puducherry, multiple R&D centres, and 17 offices across India; Schiller makes advanced healthcare equipment accessible through a network of 100 sales and service dealers across more than 45 locations. Our product range includes Critical Care, Anaesthesia, Emergency Care, Cardiology, Respiratory Diagnostics, Radiology, and Robotics.

Scan here
to know more




Part No. B9.100003  01/26 Ver. 1.2

All registered trademarks acknowledged.

SCHILLER
The Art of Diagnostics




Swiss H.Q.: SCHILLER AG, Altgasse 68, P. O. Box 1052, CH -6341 Baar, Switzerland

Indian H.O.: Schiller Healthcare India Pvt. Ltd., Advance House, Makwana Road, Andheri (East), Mumbai-59. Tel : + 91 - 22 - 61523333

 No. 15/5 & 15/6, Vazhuthavur Road, Kurumbapet, Puducherry - 605 009, India

Regional Offices: Mumbai: 022-61523333 Email: sales.west@schillerindia.com Delhi: 011-41062067 / 09312432205 Email: sales.north@schillerindia.com Kolkata: 033-23593102 / 033-23593103 Email: sales.east@schillerindia.com Chennai: 044-28232648 / 044-28311021 Email: sales.south@schillerindia.com Ahmedabad: 079- 35337949 / 079- 35337950 Email: sales.west@schillerindia.com Puducherry: 0413- 220 2680 / 0413- 229 2940 Email: sales.south@schillerindia.com

CIN : U33110MH1997PTC111307

 sales@schillerindia.com, help@schillerindia.com  www.schillerindia.com  1800 209 8998

Follow us on :  /SchillerHealthcareIndia  /schiller_healthcare_india  /schillerhealthcareindia  /schiller_healthcare_india  @schillerindia1074