


TRUSCOPE CNS

Central Monitoring System




SCHILLER

The Art of Diagnostics

TRUSCOPE CNS

- Supports nursing group, patient group
- Waveform freezing
- Large storage capacity
- Bed switching
- Supports displays of up to 4 screens
- Viewbed screen, Multibed screen, 8H dynamic Minitrends screen, ECG screen, OxyCRG screen, Big numerics screen switching
- 24h ECG summary



PORTABLE TERMINALS REMOTE VIEW

Designed for the tablets and smartphones of your choice. It goes wherever the clinician goes and allows remote monitoring of real-time patient data, securely and conveniently. It can be accessed on any device with a web browser and an internet connection.



REVIEW

REMOTE VIEW



240h

Trend data



240h

Full disclosure



1000

NIBP events



48h

OxyCRG

240h DM trends, 720 C.O. measurement results, 1000 alarm events, 8h dynamic trend, 24h ST segment, 24h ARR, 24h ECG sum, 120s freeze waveform full disclosure.

APPLICATIONS

CLINICAL ASSISTIVE APPLICATIONS


- ❖ Automatic/Manual calibration
- ❖ 24h ECG summary, ST, QT, QTc, ARR
- ❖ EWS, GCS

MULTIPLE FUNCTIONS FOR HEALTHCARE

- ❖ Provides 6 calculation functions to help clinicians
- ❖ Extendible language translation function
- ❖ Various data reports function with ST/QT (real-time), ST (review), 24h ECG summary
- ❖ Original HL7 connection



4 screens





SCHILLER

The Art of Diagnostics

About SCHILLER

SCHILLER A. G., a company based in Switzerland, was founded by Alfred Schiller in 1974. Over the last 50 years, it has evolved into a successful Global Group with 31 subsidiaries and dealers across 30 countries. Today, SCHILLER is a world-leading manufacturer and supplier of devices for cardiopulmonary diagnostics, defibrillation and patient monitoring as well as software solutions for the medical industry.

Schiller in India was founded as a joint venture in 1997, and since then has established itself as a leader in the medical technology sector. With 450+ employees, a state-of-the-art ISO 13485-certified production centre in Puducherry, multiple R&D centres, and 17 offices across India; Schiller makes advanced healthcare equipment accessible through a network of 100 sales and service dealers across more than 45 locations. Our product range includes Critical Care, Anaesthesia, Emergency Care, Cardiology, Respiratory Diagnostics, Radiology, and Robotics.

Scan here
to know more



Part no.: B2.100007 02/2026 Ver. no.: 2.1

All registered trademarks acknowledged.

SCHILLER
The Art of Diagnostics

Swiss H.Q.: SCHILLER AG, Altgasse 68, P. O. Box 1052, CH -6341 Baar, Switzerland

Indian H.O.: Schiller Healthcare India Pvt. Ltd., Advance House, Makwana Road, Andheri (East), Mumbai-59. Tel : + 91 - 22 - 61523333

No. 15/5 & 15/6, Vazhuthavur Road, Kurumbapet, Puducherry - 605 009, India

Regional Offices: Mumbai: 022-61523333 Email: sales.west@schillerindia.com Delhi: 011-41062067 / 09312432205 Email: sales.north@schillerindia.com Kolkata: 033-23593102 / 033-23593103 Email: sales.east@schillerindia.com Chennai: 044-28232648 / 044-28311021 Email: sales.south@schillerindia.com Ahmedabad: 079- 35337949 / 079- 35337950 Email: sales.west@schillerindia.com Puducherry: 0413- 220 2680 / 0413- 229 2940 Email: sales.south@schillerindia.com

CIN : U33110MH1997PTC111307

sales@schillerindia.com, help@schillerindia.com www.schillerindia.com 1800 209 8998

Follow us on : [f](#) /SchillerHealthcareIndia [i](#) /schiller_healthcare_india [in](#) /schillerhealthcareindia [e](#) /schiller_healthcare_india [y](#) @schillerindia1074