


# TRUSCOPE M10

Compact Patient Monitor



**SCHILLER**  
The Art of Diagnostics

# TRUSCOPE M10

- ❖ A patient monitor with a concise and ergonomic design, new software architecture and cutting-edge technology.

## USER-FRIENDLY DESIGN

- ❖ The ergonomic design makes it extremely convenient for users to operate
- ❖ Unit is portable with a concealed handle
- ❖ Has a highly efficient capacitive touch screen (optional) with an HD visual experience
- ❖ Operates using easy and simple gestures
- ❖ Has a fully integrated front panel without gaps and is easy to clean



The layout of the display can be changed easily



The brightness adjusts automatically based on ambient light



The fanless design reduces the risk of cross-contamination



Has 4 hours of battery life

## SIZE AND WEIGHT

- ❖ Size : 193mm X 288mm X 236mm
- ❖ Weight : < 3kg

## DISPLAY

- ❖ Type : Color TFT LCD
- ❖ Size(diagonal) : 10.4"
- ❖ Resolution : 1024 X 600 pixels

## INTERFACE

- ❖ Wired Network : RJ45 interface
- ❖ USB : 2 Sockets
- ❖ Battery Type : Rechargeable Li-ion battery

## STANDARD ACCESSORIES


- ❖ 5-lead ECG Cable
- ❖ NIBP Hose
- ❖ Adult BP Cuff (Reusable)
- ❖ SpO<sub>2</sub> Ext Cable
- ❖ SpO<sub>2</sub> Probe Adult

## OPTIONAL ACCESSORIES

- ❖ Paediatric BP Cuff (Reusable)
- ❖ Neonatal BP Cuff (Reusable)
- ❖ Neonatal SpO<sub>2</sub> Probe
- ❖ Skin Temperature Probe
- ❖ Rectal / Nasel Probe
- ❖ Neonatal Cuff (Disposable)



## RESPIRATION RATE (from the Pleth)

Pulse oximetry is the most commonly used continuous non-invasive measurement. The TRUSCOPE M10 SpO<sub>2</sub> technology can provide additional Respiration Rate facilitating early recognition of deteriorating patient conditions leading to fewer rescue interventions. These patient monitors can also reduce the consumption of disposable accessories, saving the medical staff's valuable time.




## EWS (Early Warning Scoring)

The Early Warning Scoring (EWS) in the M series monitor is a physiologic scoring system for patient assessment such as respiratory rate, heart rate, systolic blood pressure, level-of-consciousness, body temperature, etc. The EWS can detect changes in a patient's vital signs, leading to early intervention by the rapid response teams thereby preventing critical events before they happen.



## GCS (Glasgow Coma Scale)

The GCS is a neurological scale that aims to provide a reliable and objective way of recording the state of a person's consciousness for initial as well as subsequent assessment.





# SCHILLER

The Art of Diagnostics

## About SCHILLER

SCHILLER A. G., a company based in Switzerland, was founded by Alfred Schiller in 1974. Over the last 50 years, it has evolved into a successful Global Group with 31 subsidiaries and dealers across 30 countries. Today, SCHILLER is a world-leading manufacturer and supplier of devices for cardiopulmonary diagnostics, defibrillation and patient monitoring as well as software solutions for the medical industry.

Schiller in India was founded as a joint venture in 1997, and since then has established itself as a leader in the medical technology sector. With 450+ employees, a state-of-the-art ISO 13485-certified production centre in Puducherry, multiple R&D centres, and 17 offices across India; Schiller makes advanced healthcare equipment accessible through a network of 100 sales and service dealers across more than 45 locations. Our product range includes Critical Care, Anaesthesia, Emergency Care, Cardiology, Respiratory Diagnostics, Radiology, and Robotics.

Scan here  
to know more



Part no.: B2.100019 09/2024 Ver. no.: 1

All registered trademarks acknowledged.

**SCHILLER**  
The Art of Diagnostics

Swiss H.Q.: SCHILLER AG, Altgasse 68, P. O. Box 1052, CH -6341 Baar, Switzerland

Indian H.O.: Schiller Healthcare India Pvt. Ltd., Advance House, Makwana Road, Andheri (East), Mumbai-59. Tel : + 91 - 22 - 61523333

No. 15/5 & 15/6, Vazhuthavur Road, Kurumbapet, Puducherry - 605 009, India

Regional Offices: Mumbai: 022-61523333 Email: sales.west@schillerindia.com Delhi: 011-41062067 / 09312432205 Email: sales.north@schillerindia.com Kolkata: 033-23593102 / 033-23593103 Email: sales.east@schillerindia.com Chennai: 044-28232648 / 044-28311021 Email: sales.south@schillerindia.com Ahmedabad: 079- 35337949 / 079- 35337950 Email: sales.west@schillerindia.com Puducherry: 0413- 220 2680 / 0413- 229 2940 Email: sales.south@schillerindia.com

CIN : U33110MH1997PTC111307

sales@schillerindia.com, help@schillerindia.com www.schillerindia.com 1800 209 8998

Follow us on : /SchillerHealthcareIndia /schiller\_healthcare\_india /schillerhealthcareindia /schiller\_healthcare\_india @schillerindia1074