




A-mac S23 - 3 Gas Anaesthesia Machine



SCHILLER
The Art of Diagnostics



Accurate

Facilitates the effective administration of anaesthesia and promotes smooth recovery



Adaptable

Adaptable to diverse patient demographics and surgical interventions



Reliable

SCHILLER's legacy continues to thrive with over 50 years of innovation & trust


A-mac S23 - 3 Gas



Introducing **A-mac S23 - 3 Gas**, an all-encompassing anaesthesia workstation, featuring a 3 Gas system and low flow anaesthesia capabilities, designed to cater to adults, paediatrics, and infants, and integrated with a ventilator for unparalleled operational excellence.



Breathing System

- ❖ **Durable Construction** : Reliable aluminum-magnesium alloy build
- ❖ **Enhanced Efficiency** : Fast-release APL valve for improved post-surgery resuscitation
- ❖ **Patented Design** : Long-lasting, non-consumable components
- ❖ **Heightened Sensitivity** : Ultra-light ceramic one-way valve diaphragm for paediatric anaesthesia
- ❖ **Extended Operating Time** : Large 1.8 L canister for continuous use
- ❖ **Robust & Non-consumable Differential Pressure Flow Sensor**



Display

- ❖ Experience ease of use with an 8.4" TFT touch screen and rotary knob for efficient parameter settings
- ❖ Enjoy the convenience of a mute button for silent operation for 2 minutes



Monitoring System

Intuitive Monitoring


- ❖ Effortlessly track volume and pressure parameters
- ❖ Color-coded parameters for easy identification

User-Friendly Display

- ❖ Clear layout for quick navigation
- ❖ Easy operation keys for swift setting adjustments and efficient parameter confirmation

EtCO₂ & AGM Monitoring (Optional)

- ❖ Comprehensive Monitoring
- ❖ Non-invasive & User Friendly
- ❖ Integrated alarms & alerts



Alarm System

Comprehensive Alarm System

- ❖ Covers ventilation, gas, and agent alerts
- ❖ 3-level alarm system for enhanced safety and security

Alarm Log

- ❖ Easy review of past incidents for monitoring purposes



SCHILLER

The Art of Diagnostics

About SCHILLER

SCHILLER A. G., a company based in Switzerland, was founded by Alfred Schiller in 1974. Over the last 50 years, it has evolved into a successful Global Group with 31 subsidiaries and dealers across 30 countries. Today, SCHILLER is a world-leading manufacturer and supplier of devices for cardiopulmonary diagnostics, defibrillation and patient monitoring as well as software solutions for the medical industry.

Schiller in India was founded as a joint venture in 1997, and since then has established itself as a leader in the medical technology sector. With 450+ employees, a state-of-the-art ISO 13485-certified production centre in Puducherry, multiple R&D centres, and 17 offices across India; Schiller makes advanced healthcare equipment accessible through a network of 100 sales and service dealers across more than 45 locations. Our product range includes Critical Care, Anaesthesia, Emergency Care, Cardiology, Respiratory Diagnostics, Radiology, and Robotics.

Scan here
to know more:




Part No. B11.100018  5/2026 Ver. 1.3

All registered trademarks acknowledged.

SCHILLER
The Art of Diagnostics




Swiss H.Q.: SCHILLER AG, Altgasse 68, P. O. Box 1052, CH -6341 Baar, Switzerland

Indian H.O.: Schiller Healthcare India Pvt. Ltd., Advance House, Makwana Road, Andheri (East), Mumbai-59. Tel.: + 91 - 22 - 61523333

 No. 15/5 & 15/6, Vazhuthavur Road, Kurumbapet, Puducherry - 605 009, India

Regional Offices: Mumbai: 022-61523333 Email: sales.west@schillerindia.com Delhi: 011-41062067 / 09312432205 Email: sales.north@schillerindia.com Kolkata: 033-23593102 / 033-23593103 Email: sales.east@schillerindia.com Chennai: 044-28232648 / 044-28311021 Email: sales.south@schillerindia.com Ahmedabad: 079- 35337949 / 079- 35337950 Email: sales.west@schillerindia.com Puducherry: 0413- 220 2680 / 0413- 229 2940 Email: sales.south@schillerindia.com

CIN : U33110MH1997PTC111307

 sales@schillerindia.com, help@schillerindia.com  www.schillerindia.com  1800 209 8998

Follow us on :  /SchillerHealthcareIndia  /schiller_healthcare_india  /schillerhealthcareindia  /schiller_healthcare_india  @schillerindia1074