

# FRED PA-1

Simplicity as a priority



# FRED PA-1

## FEATURES

- 200J fully automatic defibrillator
- Immediate Power On
- Immediately advises to call emergency number written on the device
- Easy 3 step operation: Start -> Analysis -> Shock
- Highly efficient Lithium batteries deliver 140 shocks at maximum energy
- Automatic energy level selection, required for efficient shock delivery
- SD Card facility to store data during intervention
- Capacity of storing 500 events and 2hrs of ECG recording storage in the SD card
- Complete user guidance by voice prompts & pictograms for easy operation
- Metronome facility for accurate CPR
- LED indications for Maintenance, Replace Electrode and Replace Battery

## EASY ACCESS

Electrode holding slot

## Call

Emergency number





## No time delay

Pre-connected electrode allows faster analysis

## Safe & Ready to use

Self-test indicator with ready-to-use indication





### Saves time

Auto-switch ON once the cover is lifted.  
Time efficient device fit to save a life



### Repeat

Repeat last instruction given to user  
in case an instruction is missed out



### Warnings

LEDs indicating status of the device



### User friendly

Easy user interface with Pictograms  
give clear 1-2-3 step guide



## About SCHILLER

SCHILLER A. G., a company based in Switzerland, was founded by Alfred Schiller in 1974. Over the last 50 years, it has evolved into a successful Global Group with 31 subsidiaries and dealers across 30 countries. Today, SCHILLER is a world-leading manufacturer and supplier of devices for cardiopulmonary diagnostics, defibrillation and patient monitoring as well as software solutions for the medical industry.

Schiller in India was founded as a joint venture in 1997, and since then has established itself as a leader in the medical technology sector. With 450+ employees, a state-of-the-art ISO 13485-certified production centre in Puducherry, multiple R&D centres, and 17 offices across India; Schiller makes advanced healthcare equipment accessible through a network of 100 sales and service dealers across more than 45 locations. Our product range includes Critical Care, Anaesthesia, Emergency Care, Cardiology, Respiratory Diagnostics, Radiology, and Robotics.

Scan here  
to know more





Part No. B4.100007  2/2026 Ver. 1.1

All registered trademarks acknowledged.

**SCHILLER**  
The Art of Diagnostics

Swiss H.Q.: SCHILLER AG, Altgasse 68, P. O. Box 1052, CH -6341 Baar, Switzerland

Indian H.O.: Schiller Healthcare India Pvt. Ltd., Advance House, Makwana Road, Andheri (East), Mumbai-59. Tel : + 91 - 22 - 61523333

  : No. 15/5 & 15/6, Vazhuthavur Road, Kurumbapet, Puducherry - 605 009, India

Regional Offices: Mumbai: 022-61523333 Email: sales.west@schillerindia.com Delhi: 011-41062067 / 09312432205 Email: sales.north@schillerindia.com Kolkata: 033-23593102 / 033-23593103 Email: sales.east@schillerindia.com Chennai: 044-28232648 / 044-28311021 Email: sales.south@schillerindia.com Ahmedabad: 079-35337949 / 079-35337950 Email: sales.west2@schillerindia.com Puducherry: 0413-2202680 / 0413-2292940 Email: sales.south@schillerindia.com

CIN : U33110MH1997PTC111307

 sales@schillerindia.com, help@schillerindia.com  www.schillerindia.com  1800 209 8998

Follow us on :  /SchillerHealthcareIndia  /schiller\_healthcare\_india  /schillerhealthcareindia  /schiller\_healthcare\_india  @schillerindia1074