

neumovent.SCHILLER



Comprehensive monitoring
High-flow oxygen therapy
New ventilation modes



Mechanical ventilator that provides
invasive and non-invasive ventilation,
for neonatal & paediatric patients



SCHILLER
The Art of Diagnostics

GraphNet neo

GraphNet neo is tailored for neonatal intensive care, offering both invasive and non-invasive support for infants, including premature babies, while closely monitoring key respiratory parameters.




COMPREHENSIVE MONITORING

- The proximal flow sensor enhances flow and volume monitoring in the patient circuit, enabling accurate ventilation parameter measurements for informed decision-making, effective treatment, and ensured patient safety



Because they deserve the best care.






- ❖ Ventilation for neonatal patients
- ❖ Available for both invasive and non-invasive ventilation
- ❖ Continuous flow CPAP: non-invasive ventilation using nasal prongs with automatic leak compensation
- ❖ High-flow oxygen therapy
- ❖ Proximal flow sensor for enhanced monitoring
- ❖ Respiratory mechanics menu
- ❖ Suitable for intra-hospital transportation
- ❖ Accurate tidal volumes starting from 2mL
- ❖ 72 hours of trend storage
- ❖ Built-in battery with over 2.5 hours capacity
- ❖ Numerical display screen

TARGETED VENTILATION MODE

- ❖ PRVC is an adaptive mode that automatically adjusts the pressure setpoint to ensure the user-configured volume, while accounting for potential changes in the patient's lung impedance



HIGH-FLOW OXYGEN THERAPY

- ❖ The addition of this method expands the use of the ventilator to the stages that are prior to and after the mechanical ventilation itself, by providing non-invasive assistance to oxygenation




About SCHILLER

SCHILLER A. G., a company based in Switzerland, was founded by Alfred Schiller in 1974. Over the last 50 years, it has evolved into a successful Global Group with 31 subsidiaries and dealers across 30 countries. Today, SCHILLER is a world-leading manufacturer and supplier of devices for cardiopulmonary diagnostics, defibrillation and patient monitoring as well as software solutions for the medical industry.

Schiller in India was founded as a joint venture in 1997, and since then has established itself as a leader in the medical technology sector. With 450+ employees, a state-of-the-art ISO 13485-certified production centre in Puducherry, multiple R&D centres, and 17 offices across India; Schiller makes advanced healthcare equipment accessible through a network of 100 sales and service dealers across more than 45 locations. Our product range includes Critical Care, Anaesthesia, Emergency Care, Cardiology, Respiratory Diagnostics, Radiology, and Robotics.

Scan here
to know more:





Part no.: B11.100015  12/24 Ver. no.: 1

All registered trademarks acknowledged. Accessories shown are not part of standard configuration

SCHILLER
The Art of Diagnostics




Swiss H.Q.: SCHILLER AG, Altgasse 68, P. O. Box 1052, CH -6341 Baar, Switzerland

Indian H.O.: Schiller Healthcare India Pvt. Ltd., Advance House, Makwana Road, Andheri (East), Mumbai-59. Tel: + 91 - 22 - 61523333

  : No. 15/5 & 15/6, Vazhuthavur Road, Kurumbapet, Puducherry - 605 009, India

Regional Offices: Mumbai: 022-61523333 Email: sales.west@schillerindia.com Delhi: 011-41062067 / 09312432205 Email: sales.north@schillerindia.com Kolkata: 033-23593102 / 033-23593103 Email: sales.east@schillerindia.com Chennai: 044-28232648 / 044-28311021 Email: sales.south@schillerindia.com Ahmedabad: 079-35337949 / 079-35337950 Email: sales.west2@schillerindia.com Puducherry: 0413-2202680 / 0413-2292940 Email: sales.south@schillerindia.com

CIN : U33110MH1997PTC111307

 sales@schillerindia.com, help@schillerindia.com  www.schillerindia.com  1800 209 8998

Follow us on :  /SchillerHealthcareIndia  /schiller_healthcare_india  /schillerhealthcareindia  /schiller_healthcare_india  @schillerindia1074