

MTM-1500 MED MTM-1400 MED

High-quality, low-maintenance treadmills



SCHILLER
The Art of Diagnostics

SCHILLER
MTM

MTM-1500 MED

SAFE – ACCURATE – QUIET


SCHILLER's treadmill ergometer MTM-1500 med is designed for exercise tests in sports medicine, cardiology, physiotherapy and fitness training. The ergometer is equipped and tested with the latest safety systems and has been manufactured to meet stringent quality requirements.

KEY FEATURES:

- ❖ Compatible with all SCHILLER exercise ECG devices
- ❖ With potential equalisation transformer for potential isolation from the mains, according to IEC 60601-1
- ❖ Automatic control by computer via RS-232 interface (interface cable 5 m included)
- ❖ Individually adjustable emergency stop button
- ❖ Steel tube hand rail; can be mounted on the right or left hand side. A second hand rail is optionally available.
- ❖ Automatic centring mechanism for the belt

ADVANTAGES:

- ❖ Very high durability thanks to outstanding manufacturing quality
- ❖ Excellent price / performance ratio
- ❖ Extremely safe
- ❖ Almost noiseless drive
- ❖ Maintenance-free motor
- ❖ Modern design
- ❖ Generation of reports with terminal software



MTM-1400 MED

THE NEW SCHILLER TREADMILL MTM-1400 MED

STABLE AND LOW-MAINTENANCE

The MTM-1400 med treadmill is virtually indestructible with its sturdy frame. It is very low maintenance and provides the runner or patient a comfortable run with its advanced construction. MTM-1400 med is distinguished by its smooth running, various functions, strong drive engine and the timeless design.

MEDICAL DEVICE TREADMILL

As a medical device, it is particularly suitable for use in the fields of cardiology, neurology, cardiac rehabilitation and physiotherapy.

The interface of the SCHILLER MTM-1400 med treadmill ergometer allows connection to ECG, ergospirometry systems, blood pressure measuring devices and software programs via the coscom® protocol.

CHILDREN HANDRAILS FOR PAEDIATRIC REHABILITATION

Children handrails are available as an additional option (only in combination with the safety arch). The MTM-1400 med treadmill can be used for paediatric rehabilitation with this option. Children can hold onto the handrail both on the sides and in front while walking on the treadmill, and therefore always have a sense of security and control.

NUMEROUS ADDITIONAL OPTIONS

You can customise the MTM-1400 med treadmill to your field of application – whether rehabilitation, cardiopulmonary diagnostics for patients or fitness – with the numerous additional options.





About SCHILLER

SCHILLER A. G., a company based in Switzerland, was founded by Alfred Schiller in 1974. Over the last 50 years, it has evolved into a successful Global Group with 31 subsidiaries and dealers across 30 countries. Today, SCHILLER is a world-leading manufacturer and supplier of devices for cardiopulmonary diagnostics, defibrillation and patient monitoring as well as software solutions for the medical industry.

Schiller in India was founded as a joint venture in 1997, and since then has established itself as a leader in the medical technology sector. With 450+ employees, a state-of-the-art ISO 13485-certified production centre in Puducherry, multiple R&D centres, and 17 offices across India; Schiller makes advanced healthcare equipment accessible through a network of 100 sales and service dealers across more than 45 locations. Our product range includes Critical Care, Anaesthesia, Emergency Care, Cardiology, Respiratory Diagnostics, Radiology, and Robotics.

Scan here
to know more:




Part No. B9.100024  01/2026 Ver no . 1.2

All registered trademarks acknowledged.

SCHILLER
The Art of Diagnostics




Swiss H.Q.: SCHILLER AG, Altgasse 68, P. O. Box 1052, CH -6341 Baar, Switzerland

Indian H.O.: Schiller Healthcare India Pvt. Ltd., Advance House, Makwana Road, Andheri (East), Mumbai-59. Tel : + 91 - 22 - 61523333

 : No. 15/5 & 15/6, Vazhuthavur Road, Kurumbapet, Puducherry - 605 009, India

Regional Offices: Mumbai: 022-61523333 Email: sales.west@schillerindia.com Delhi: 011-41062067 / 09312432205 Email: sales.north@schillerindia.com Kolkata: 033-23593102 / 033-23593103 Email: sales.east@schillerindia.com Chennai: 044-28232648 / 044-28311021 Email: sales.south@schillerindia.com Ahmedabad: 079-35337949 / 079-35337950 Email: sales.west2@schillerindia.com Puducherry: 0413-2202680 / 0413-2292940 Email: sales.south@schillerindia.com

CIN : U33110MH1997PTC111307

 sales@schillerindia.com, help@schillerindia.com  www.schillerindia.com  1800 209 8998

Follow us on :  /SchillerHealthcareIndia  /schiller_healthcare_india  /schillerhealthcareindia  /schiller_healthcare_india  @schillerindia1074