


TRUSCOPE M10

Experience the future of patient monitoring with a sleek, ergonomic design and advanced technology. It features an optional high-definition capacitive touchscreen for a seamless, intuitive user experience.



Standard Parameters

- ❖ ECG
- ❖ SpO₂
- ❖ NIBP
- ❖ RR
- ❖ Temp

- ❖ **Advanced Monitoring:** ST Analysis, ST Mapping, QT analysis, 27 types of arrhythmia analysis, GCS & EWS
- ❖ **Comprehensive Calculations:** Drug dose, hemodynamic, renal, oxygenation and ventilation profiles
- ❖ **OXY CRG Software:** For detailed respiratory monitoring



The layout of the display can be changed easily



The fanless design reduces the risk of cross-contamination



The brightness adjusts automatically based on ambient light




Has 4 hours of battery life

DISPLAY

- ❖ Type : Color TFT LCD
- ❖ Screen Size(diagonal) : 10.4"
- ❖ Resolution : 1024 X 600 pixels
- ❖ Weight : < 3kg

INTERFACE

- ❖ Wired Network : RJ45 interface
- ❖ USB : 2 Sockets
- ❖ Battery Type : Rechargeable Li-ion battery
- ❖ Touchscreen (Optional)



About SCHILLER

SCHILLER A. G., a company based in Switzerland, was founded by Alfred Schiller in 1974. Over the last 50 years, it has evolved into a successful Global Group with 31 subsidiaries and dealers across 30 countries. Today, SCHILLER is a world-leading manufacturer and supplier of devices for cardiopulmonary diagnostics, defibrillation and patient monitoring as well as software solutions for the medical industry.

Schiller in India was founded as a joint venture in 1997, and since then has established itself as a leader in the medical technology sector. With 450+ employees, a state-of-the-art ISO 13485-certified production centre in Puducherry, multiple R&D centres, and 17 offices across India; Schiller makes advanced healthcare equipment accessible through a network of 100 sales and service dealers across more than 45 locations. Our product range includes Critical Care, Anaesthesia, Emergency Care, Cardiology, Respiratory Diagnostics, Radiology, and Robotics.

Scan here
to know more



Part No.: B2.100019 07/2026 Ver. no.: 1.2

All registered trademarks acknowledged.

SCHILLER
The Art of Diagnostics

Swiss H.Q.: SCHILLER AG, Altgasse 68, P. O. Box 1052, CH -6341 Baar, Switzerland

Indian H.O.: Schiller Healthcare India Pvt. Ltd., Advance House, Makwana Road, Andheri (East), Mumbai-59. Tel : + 91 - 22 - 61523333

Plot : No. 15/5 & 15/6, Vazhuthavur Road, Kurumbapet, Puducherry - 605 009, India

Regional Offices: Mumbai: 022-61523333 Email: sales.west@schillerindia.com Delhi: 011-41062067 / 09312432205 Email: sales.north@schillerindia.com Kolkata: 033-23593102 / 033-23593103 Email: sales.east@schillerindia.com Chennai: 044-28232648 / 044-28311021 Email: sales.south@schillerindia.com Ahmedabad: 079-35337949 / 079-35337950 Email: sales.west2@schillerindia.com Puducherry: 0413-2202680 / 0413-2292940 Email: sales.south@schillerindia.com

CIN : U33110MH1997PTC111307

sales@schillerindia.com, help@schillerindia.com www.schillerindia.com 1800 209 8998

Follow us on : /SchillerHealthcareIndia /schiller_healthcare_india /schillerhealthcareindia /schiller_healthcare_india @schillerindia1074