

DEFIGARD 400

Portable Biphasic Defibrillator



SCHILLER

The Art of Diagnostics

DEFIGARD 400

DEFIGARD 400 is a versatile biphasic defibrillator. Specifically designed for critical care, it combines portability, precision and power.

- ✧ A compact biphasic defibrillator designed for a quick and efficient emergency response for both adult and paediatric patients.
- ✧ With its advanced features and intuitive design, it's built for optimal performance in critical care situations.



ENHANCED USABILITY AND PERFORMANCE



Large Display: 8.4" TFT color display ensures clear visibility in any condition



User-Friendly Interface: Rotary knob for quick settings adjustments and easy navigation



Versatile Application: Suitable for both adult and paediatric patients and is ideal for various clinical environments



Memory: Stores the last 11 defibrillation events and 24 hours of detailed parameter trends for effective monitoring and analysis with option of additional USB storage



Extended Monitoring: Capable of 24 hours of tabular trend monitoring for all parameters



Performance: Supports 2-3 hours of continuous monitoring and delivers up to 100 shocks at 200J energy level

Sync Mode: Synchronized defibrillation mode to deliver the shock at the precise point in the cardiac cycle ensuring restoration of normal sinus rhythm without complications.


Contact Indicator: Paddles contact indicator present on screen for good paddle contact. If Paddles are not connected to the patient properly, then the poor contact identification is displayed.

Safety Discharge: If the electrodes are not properly attached to the patient, pressing the 'Discharge' button won't shock the patient. Instead, the defibrillator will internally discharge the energy. Also, if the defibrillator stays charged for more than 20 seconds without being used, it will safely discharge energy internally.

EXTENDED FUNCTIONALITIES FOR VERSATILE USE

Enhance the DEFIGARD 400 with optional upgrades to tailor the device to specific clinical needs:

- ✧ **Automated External Defibrillator (AED) Mode:** Provides guided defibrillation with visual instructions for easy use by responders of all skill levels.
- ✧ **Integrated Thermal Printer:** Offers both automatic and manual printing for precise intervention records.
- ✧ **USB Interface for Data Storage:** Enables efficient patient data storage, streamlining workflow and data management.
- ✧ **Internal Defibrillation:** Internal defibrillation mode for open-heart surgical procedures, with a full range of spoon sizes and energy levels.
- ✧ **Advanced Monitoring Capabilities:** 360 J, EtCO₂ (Mainstream), real-time CPR feedback, external pacemaker, 12-lead ECG, SpO₂ (Digital / Nellcor) and NIBP for thorough patient evaluation.



About SCHILLER

SCHILLER A. G., a company based in Switzerland, was founded by Alfred Schiller in 1974. Over the last 50 years, it has evolved into a successful Global Group with 31 subsidiaries and dealers across 30 countries. Today, SCHILLER is a world-leading manufacturer and supplier of devices for cardiopulmonary diagnostics, defibrillation and patient monitoring as well as software solutions for the medical industry.

Schiller in India was founded as a joint venture in 1997, and since then has established itself as a leader in the medical technology sector. With 450+ employees, a state-of-the-art ISO 13485-certified production centre in Puducherry, multiple R&D centres, and 17 offices across India; Schiller makes advanced healthcare equipment accessible through a network of 100 sales and service dealers across more than 45 locations. Our product range includes Critical Care, Anaesthesia, Emergency Care, Cardiology, Respiratory Diagnostics, Radiology, and Robotics.

Scan here
to know more:



Part no.: B1.100007  07/2026 Ver. no.: 1.5

All registered trademarks acknowledged.

SCHILLER
The Art of Diagnostics

Swiss H.Q.: SCHILLER AG, Altgasse 68, P. O. Box 1052, CH - 6341 Baar, Switzerland

Indian H.O.: Schiller Healthcare India Pvt. Ltd., Advance House, Makwana Road, Andheri (East), Mumbai-59. Tel : + 91 - 22 - 61523333

 : No. 15/5 & 15/6, Vazhuthavur Road, Kurumbapet, Puducherry - 605 009, India

Regional Offices: Mumbai: 022-61523333 Email: sales.west@schillerindia.com Delhi: 011-41062067 / 09312432205 Email: sales.north@schillerindia.com Kolkata: 033-23593102 / 033-23593103 Email: sales.east@schillerindia.com Chennai: 044-28232648 / 044-28311021 Email: sales.south@schillerindia.com Ahmedabad: 079-35337949 / 079-35337950 Email: sales.west2@schillerindia.com Puducherry: 0413-220 2680 / 0413-229 2940 Email: sales.south@schillerindia.com

CIN : U33110MH1997PTC111307

 sales@schillerindia.com, help@schillerindia.com  www.schillerindia.com  1800 209 8998

Follow us on :  /SchillerHealthcareIndia  /schiller_healthcare_india  /schillerhealthcareindia  /schiller_healthcare_india  @schillerindia1074