

TRUSCOPE P1


Portable Transport Monitor



SCHILLER
The Art of Diagnostics

TRUSCOPE P1 - Portable Transport Monitor

The TRUSCOPE P1 transport monitor ensures continuous patient monitoring throughout every stage of care, enabling a seamless transition from ambulance to operating room, ICU, inter-intra hospital transport and recovery. Its compact, portable design helps maintain uninterrupted access to vital patient data wherever care is desired.



Docking Station

Transfer Handle

- ❖ Lightweight and compact (< 0.9 kg)
- ❖ Color TFT 5.5" HD Capacitive touchscreen
- ❖ Long battery life over 8 hours
- ❖ Drop resistant
- ❖ Dust & liquid ingress protection
- ❖ Supports gesture control for various screen option & controls

- ❖ **Advanced Monitoring:** ST Analysis, ST Mapping, QT & 27 types of arrhythmia analysis.
- ❖ **Comprehensive Calculations:** Drug dose, hemodynamic, renal, oxygenation and ventilation profiles.
- ❖ **OXY CRG Software:** For detailed respiratory monitoring.
- ❖ **Transport Options:** Transfer handle, docking station, transfer bracket (for ambulance) & AC DC Adapter.
- ❖ Bed rail option to hook monitor to Bed rails etc.

Adaptive Monitoring Across Every Care Setting

Seamless patient monitoring with modular flexibility—designed for both indoor care settings and outdoor emergencies.



Operation Room

Ensures continuous and accurate vital monitoring during surgical procedures, supporting clinical decision-making throughout the operation.

Recovery Room

Supports safe recovery with stable monitoring of vital parameters during post-surgical observation.

ICU

Advanced multi-parameter monitoring ensures continuous, high dependency patient assessment with accurate data.

Outdoor Emergency

Compact and lightweight module design allows quick deployment for immediate patient assessment in critical situations.

Ambulance

Ensures uninterrupted monitoring with long battery backup and secure data continuity during intra- and inter-hospital transfers.



TRUSCOPE P1 as a Multiparameter Module

- ✦ Slides seamlessly into the module slot of the P-Series Monitor, enabling an expanded range of parameter measurements
- ✦ Connects to the Central Monitoring System over a wireless network for an untethered and seamless monitoring experience



Touch-enabled Portable device
5.5 inch



Easy to Carry
<0.9kg



Efficient Power
8 hours



Additive module with Docking station
Masimo/Nellcor, SpO₂ Sidestream /
Mainstream CO₂, BIS, EEG, NMT, ICG,
Dual IBP or AG



Medical Transport cket




Bed Rail clamp for Docking Station / Transfer Handle

Available Configurations



- P1 Standard configuration: 3/5/6-lead ECG, RESP, NIBP, TEMP, SpO₂
- P1 Standard configuration+12-lead
- P1 Standard configuration+IBP
- P1 Standard configuration+IBP+12-lead



About SCHILLER

SCHILLER A. G., a company based in Switzerland, was founded by Alfred Schiller in 1974. Over the last 50 years, it has evolved into a successful Global Group with 31 subsidiaries and dealers across 30 countries. Today, SCHILLER is a world-leading manufacturer and supplier of devices for cardiopulmonary diagnostics, defibrillation and patient monitoring as well as software solutions for the medical industry.

Schiller in India was founded as a joint venture in 1997, and since then has established itself as a leader in the medical technology sector. With 450+ employees, a state-of-the-art ISO 13485-certified production centre in Puducherry, multiple R&D centres, and 17 offices across India; Schiller makes advanced healthcare equipment accessible through a network of 100 sales and service dealers across more than 45 locations. Our product range includes Critical Care, Anaesthesia, Emergency Care, Cardiology, Respiratory Diagnostics, Radiology, and Robotics.

Scan here
to know more



All registered trademarks acknowledged.

Part no.: B2.100020 08/2026 Ver. no.: 1

SCHILLER
The Art of Diagnostics

Swiss H.Q.: SCHILLER AG, Altgasse 68, P. O. Box 1052, CH -6341 Baar, Switzerland

Indian H.O.: Schiller Healthcare India Pvt. Ltd., Advance House, Makwana Road, Andheri (East), Mumbai-59. Tel : + 91 - 22 - 61523333

India : No. 15/5 & 15/6, Vazhuthavur Road, Kurumbapet, Puducherry - 605 009, India

Regional Offices: Mumbai: 022-61523333 Email: sales.west@schillerindia.com Delhi: 011-41062067 / 09312432205 Email: sales.north@schillerindia.com Kolkata: 033-23593102 / 033-23593103 Email: sales.east@schillerindia.com Chennai: 044-28232648 / 044-28311021 Email: sales.south@schillerindia.com Ahmedabad: 079- 35337949 / 079- 35337950 Email: sales.west2@schillerindia.com Puducherry: 0413- 220 2680 / 0413- 229 2940 Email: sales.south@schillerindia.com

CIN : U33110MH1997PTC111307

sales@schillerindia.com, help@schillerindia.com www.schillerindia.com 1800 209 8998

Follow us on : /SchillerHealthcareIndia /schiller_healthcare_india /schillerhealthcareindia /schiller_healthcare_india @schillerindia1074